

## Ligation Protocol with T4 DNA Ligase (M0202)

Protocols.io also provides an [interactive version of this protocol](#) where you can discover and share optimizations with the research community.

### Protocol

1. Set up the following reaction in a microcentrifuge tube on ice.  
(*T4 DNA Ligase should be added last. Note that the table shows a ligation using a molar ratio of 1:3 vector to insert for the indicated DNA sizes.*) Use [NEBioCalculator](#) to calculate molar ratios.

COMPONENT	20 µl REACTION
10XT4 DNA Ligase Buffer*	2 µl
Vector DNA (4 kb)	50 ng (0.020 pmol)
Insert DNA (1 kb)	37.5 ng (0.060 pmol)
Nuclease-free water	to 20 µl
T4 DNA Ligase	1 µl

\* *The T4 DNA Ligase Buffer should be thawed and resuspended at room temperature.*

2. Gently mix the reaction by pipetting up and down and microfuge briefly.
3. For cohesive (sticky) ends, incubate at 16°C overnight or room temperature for 10 minutes.
4. For blunt ends or single base overhangs, incubate at 16°C overnight or room temperature for 2 hours (*alternatively, high concentration T4 DNA Ligase can be used in a 10 minute ligation*).
5. Heat inactivate at 65°C for 10 minutes.
6. Chill on ice and transform 1-5 µl of the reaction into 50 µl competent cells.